

Sample Type

➤ Types of sample "Sputum"

➤ Cryo box container holding 50ml falcon tube

➤ Two large transportation boxes for specimen with thermometer.




Referral Lab for sample s

Districts	Referral Lab
Panchthar, IlamJhapa, Morang Sunsri, Tehrathum, Sakhuwasabha, Saptari, Udayapur, Siraha,Dhanusha, Mahottari, SarlihiRautahat, Bara	IOM , Damak, Jhapa
Kathmanadu, Bhaktapur, Lalitpur,S.Palchok. Dolakha, Ramechhap,S.KhambuKavre, Nowakot, Dhading, Tanahun,Syanja ,Baglung, Kaski,Chitawan, MakawanpurParsa	NATA, Kathmandu
Darchula, Dadeldhura, Achham, Bajura,Doti,Kailali, Kanchanpur, Bardiya,Banke, Jumla, Dailakha, Surkhet, Rukum, Salyan, Dang, Kapilvastu, Arghakhasnchi, Gulmi, Palpa, Rupandehi, Nawalparasi	NTC, Bhaktapur



Sample Transportation Mechanism

- ❖ During the field operation messenger will be responsible for transporting the sputum sample to referral Lab.
- ❖ Samples will be transported with in 2 days from the camp sites.



Major support of RTLO/DTLO

- During the field operation provide support in order to maintain cold chain for samples.
 - Coordinate hospitals/local health facilities
 - Provide required support during sample transportation.



Remarks/feedback

